

DIMAR

Medical Devices

DimAir Line

Intensive Care



CPAP
PAEDIATRIC and NEONATAL HELMET

Single Patient use Helmets for the treatment of respiratory failure in CPAP and BI-PAP techniques in Neonatal Patients from 3 to 5 Kg and Paediatric Patients from 5 to 15 Kg. Supplied with a pneumatic sealed zip, they allow a quick and easy access to the patient. Comfortable and practical.

TECHNICAL FEATURES

- Opening zip with pneumatic seal
- Made of biocompatible plastic material
- 4 symmetrical Right and Left connectors for a practical and easy connection with the equipment placed on any side
- Soft Pad fabricated from medical approved anti-bacterial CoolMax® material
- Sealed access ports from 4 to 19 mm
- 2 fixing eyelets on the top of the helmet for an additional connection
- 22M, 22F ISO connectors
- Single Patient use
- Latex Free



The additional "eyelet" system allows to fix the helmet to the bed

ADVANTAGES

- Quick access to Patient's face
- Possibility to carry out cleaning, feeding and weaning operations without taking the helmet off
- Extremely Comfortable, Practical and Effective
- It avoids the infections caused by oral-nasal tracheal intubation
- Devoid of rigid ring, it prevents the formation of ulcerations to the Patient's neck
- Gases inlet and outlet can be inverted or connected on the same side; in this way Patients can also be placed on their side
- It avoids necrosis or ulcerations caused by facial masks or nasal pillows
- Associated to PEEP DIMAR Valve, it assures the two-way flow in case of lack of pressure, preventing asphyxia risk

REFERENCES

MODEL	CODE	SIZE	Internal Volume (L)
Paediatric	500/5448	Ø 60 mm	7,5 Litres
Neonatal	500/5447	Ø 50 mm	6,2 Litres

- Single pack
- Multiple pack with 5 pieces